

Clustering Sheet for Flowcell

AAFKVYYM5

<b>Name</b> Bartlomiej Gebarski	<b>Sequencer</b> NextSeq2k1	<b>Run Type</b> NextSeq2000 P2 SR100 (r1: 100, r2: 0)
<b>Clustering Date</b> 2024-05-22	<b>Sequencing Date</b> 2024-05-22	<b>Status</b> Sequenced
<b>1st Index Length</b> 15	<b>2nd Index Length</b> 15	
<b>Comments</b> no comment	<b>Reagents Checked Out</b> <input type="checkbox"/>	

- NextSeq 1000/2000 P2 Reagents (100 Cycles) v3: 1

num	Id	Scientist	Custom	PR	Comments
1	M18572_R17198_1	<div></div>	<input type="checkbox"/>	Standard	GLOE-seq

## Instructions

1. Check, if instrument is idle. If run just finished - remove used cartridge (follow instructions on the screen).
2. Prepare Cluster Sheet, Sample and all reagents.
3. Take out the flowcell from the fridge to equilibrate it to room temperature (5-10min).
4. Install Flowcell in the cartridge
5. Mix all reagents in the Eppendorf tube:
  - **1.25µl** of **1nM** PhiX;
  - **9.5µl** of **2.5nM** Library;
  - **14.25µl** of Dilution Buffer;
6. Vortex and spin down.
7. Transfer **20µl** into NextSeq 2000 cartridge in well marked 'sample'.
8. Load reagents into instrument, following instructions on the screen.